

Transscleral Cyclophotocoagulation Probe

Footplate designed for precise placement around the circumference of the limbus



Remove footplate to use MicroPulse® settings*

Part Number	Laser Compatibility
D4700A	Iridex iQ 810™, Oculight® SL, Oculight® SLx

- ◀ Designed to work with an infrared laser
- ◀ Designed to work for two different procedures
- ◀ Sold in boxes of 5

Contact your local rep to schedule a surgical evaluation or call 888-4KATALYST(888-452-8259) or 636-536-5950 to place an order!

*Iridex® and MicroPulse® are registered US trademarks of Iridex Corporation